

Dialysis or diafiltration?

D-Tube™ dialyzers and Amicon® Ultra filters give YOU the choice.

Each protein preparation is unique. Give it the special treatment it deserves with a perfectly designed device for buffer exchange, desalting, and removing solutes like urea and detergents. Select between fast and gentle dialysis using the D-Tube™ Dialyzers or diafiltration using the faster Amicon® Ultra Centrifugal Filters.

What's the difference?

Dialysis is the traditional method for desalting (removing microsolute) or buffer/solvent exchange, using osmotic pressure to drive solutes across a membrane.

Diafiltration, or ultrafiltration, achieves desalting or buffer exchange through the use of centrifugal force or other external pressure to drive small microsolute through a porous membrane. The membrane does not allow macrosolute (larger than the pore size) to pass through.



D-Tube™ dialyzers enable you to:

- Prevent precipitation or over-concentration of intractable or sensitive samples.
- Complete dialysis in just two to five hours, thanks to the D-Tube™ dialyzer's double membrane, which spreads the sample over a large surface area.
- Easily handle samples in a capped centrifuge tube format with dialysis membrane windows.
- Achieve >97% sample volume recovery from solution.

Amicon® Ultra centrifugal filters enable you to:

- Rapidly desalt or exchange buffers without diluting your sample – just concentrate the sample, discard filtrate, and reconstitute the concentrate with desired solvent.
- Repeat “washing out” until solvent exchange or desalting is complete.
- Efficiently achieve salt (microsolute) transfer across the membrane, regardless of microsolute concentration or size as quick as 30 minutes.

Ask your sample what it needs.

Based on your sample preparation needs, use the chart below to determine whether dialysis or diafiltration is the best solution.

How to choose D-Tubes™ dialyzers vs. Amicon® Ultra filters for dialysis and diafiltration

Need	D-Tube™	Amicon® Ultra
Faster access to sample		X
Sensitive samples or samples which may precipitate at higher concentrations	X	
Robust samples able to be concentrated		X
No concentration needed	X	
Concentration needed		X
Limited amounts of exchange solvent		X
Sample needs to be kept cold	Cold temperature reduces speed	Minimal effect of cold temperature on speed

Choose an Amicon® Ultra filter:

Simply use the tables below to first identify the best Molecular Weight Cut-Off (MWCO) for your sample and then select the Amicon® filter that best suits your application based on starting and ideal final volume.

Proteins	Nanoparticles	Single- and Double-Stranded Nucleic Acids	MWCO
6 < MW < 20 k	1.5 nm < dia < 3 nm		3,000
20 < MW < 60 k	3 nm < dia < 5 nm		10,000
60 < MW < 100 k	5 nm < dia < 7 nm	137-1159 bp	30,000
100 < MW < 200 k	7 nm < dia < 10 nm		50,000
200 k < MW	10 nm < dia		100,000

Amicon® Ultra Centrifugal Filters

Product	Amicon® Ultra-0.5	Amicon® Ultra-2	Amicon® Ultra-4	Amicon® Ultra-15
Maximum initial sample volume (mL)	0.5	2	4	15
Final concentrate (retentate) volume (µL)	15-20	15-70	30-70	150-300
MWCO	Qty/Pk			
3,000	8	UFC500308		UFC800308
MWCO	24	UFC500324	UFC200324	UFC800324
	96	UFC500396		UFC800396
	500	UFC5003BK		UFC900396
10,000	8	UFC501008		UFC801008
MWCO	24	UFC501024	UFC201024	UFC801024
	96	UFC501096		UFC801096
	500	UFC5010BK		UFC901096
30,000	8	UFC503008		UFC803008
MWCO	24	UFC503024	UFC203024	UFC803024
	96	UFC503096		UFC803096
	500	UFC5030BK		UFC903096
50,000	8	UFC505008		UFC805008
MWCO	24	UFC505024	UFC205024	UFC805024
	96	UFC505096		UFC805096
	500	UFC5050BK		UFC905096
100,000	8	UFC510008		UFC810008
MWCO	24	UFC510024	UFC210024	UFC810024
	96	UFC510096		UFC810096
	500	UFC5100BK		UFC910096

Amicon® Ultra filter selector tool

Find the right filter to concentrate your sample – search with our online Amicon® Ultra Selector Tool.

Visit: www.millipore.com/FastEasy



www.merckmillipore.com/offices

Merck Millipore and the M mark are trademarks of Merck KGaA, Darmstadt, Germany. D-Tube is a trademark of Merck KGaA, Darmstadt, Germany. Amicon is a registered trademark of Merck KGaA, Darmstadt, Germany.
Lit No. PR3888ENEU LS SBU-12-06314 4/2012 Printed in the USA.
© 2012 EMD Millipore Corporation, Billerica, MA USA. All rights reserved.

Choose a D-Tube™ dialyzer:

What is the starting volume of your solution? What is the molecular weight of your molecule of interest? Answer these questions, and you're ready to choose a D-Tube™ Dialyzer with the right molecular weight cut-off (MWCO)

Product	D-Tube™ Mini	D-Tube™ Midi	D-Tube™ Maxi	D-Tube™ Mega	D-Tube™ Mega	D-Tube™ Mega
Maximum initial sample volume	10 to 250 µL	50 to 800 µL	100 µL to 3 mL	3 to 10 mL	10 to 15 mL	15 to 20 mL
MWCO	Qty/Pk					
3,500	10		71506-3	71508-3	71739-3	71742-3
	50				71739-4	71742-4
6,000-8,000	10	71504-3	71507-3	71509-3	71740-3	71743-3
	50				71740-4	71743-4
	1 plate of 96	71712-3				71746-4
12,000-14,000	10	71505-3		71510-3		
	1 plate of 96	71713-3				
Floating Rack						
	10	71512-3	71513-3	71514-3	71748-3	71748-3

To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 01805 045 645

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: www.merckmillipore.com/offices

For Technical Service visit:

www.merckmillipore.com/techservice

Get Connected!

Join Merck Millipore Bioscience on your favorite social media outlet for the latest updates, news, products, innovations, and contests!

facebook.com/MerckMilliporeBioscience

twitter.com/Merck4Bio